

Buffer Preparation (Gozani Lab)

1. 1 M Tris-HCl Buffers

pH	Volume (L)	TrisBase (g)	HCl (ml)
pH 7.0	2	242.2	150-155
pH 7.5	2	242.2	120-125
pH 8.0	2	242.2	80-85

Autoclavable.

2. EDTA 0.5 M (pH8.0)

0.5M, 1L: 148 g EDTA

+ ~30-40 g NaOH to adjust pH

(or 186 g EDTA-Na.2H₂O + ~20 g NaOH)

Note: pH adjusted by NaOH is essential for solubility. Autoclavable.

3. TAE DNA Electrophoresis Buffer (50 X)

(2 M Tris, 50 mM EDTA)

2 L

484 g Tris

114.2 ml glacial acetic acid

200 ml 0.5 M EDTA 8.0

To make 1x TAE 20 L, add 400 ml 50X buffer into 19.6 L ddH₂O.

4. SDS-PAGE Gel Solutions

	Vol (L)	Tris (g)	HCl (ml)	10% SDS (ml)
4x Lower gel buffer <i>1.5 M Tris-Cl, pH 8.8, 0.4% SDS</i>	2	363.3	50-60	80 ml
4x Upper gel buffer <i>0.5 M Tris-Cl, pH 6.8, 0.4% SDS</i>	2	121.1	70-80	80 ml

4.1 10% SDS

1L:

100g SDS into 1 L, heat to 68°C for solubility. pH ~6.6.

5. 5X SDS Loading Sample Buffer

100 ml

	Stock solution	Add volume
250 mM TrisHCl pH6.8	2 M	12.5 ml
10% SDS		10 g
30% Glycerol		30 ml
5% β -mercapitoethanol (or 0.5M DTT)		5 ml
0.02% bromophenol blue	0.04%	52 ml

6. 6X DNA loading sample buffer:

(40% sucrose, 0.01-0.02% BPB)

100 ml

Add 40 g sucrose to 50 ml 0.04% BPB solution, adjust final volume 100 ml.

7. SDS-PAGE Electrophoresis Running Buffer (10x)

(1x: 25 mM Tris, 192 mM glycine, 0.1% SDS, pH8.3)

10 L.

303 g Trisbase (FW 121.1)

1440 g glycine (FW 75.07)

100 g SDS

No need to adjust pH

8. Transfer Buffer without SDS (10x)

(1x: 25 mM Tris, 192 mM glycine, pH8.3)

10 L

303 g Trisbase,

1440 g glycine

No need to adjust pH

8.1 Transfer Buffer (1x)

500 ml

50 ml of 10x Transfer buffer (without SDS) or 10x SDS-PAGE running buffer (w/ SDS)

100 ml of Methanol (final **20% methanol**)

350 ml ddH₂O

9. TBS (10x)

(1x: 150 mM NaCl, 10 mM Tris pH8.0)

10 L
876.6 g NaCl (FW 58.44),
121.1 g Tris,
~40 ml HCl
to pH8.0

9.1 TBS-T (1x)

1L
100 ml 10x TBS
10 ml 10% Tween20 (final 0.1% v/v)
890 ml ddH₂O

9.2 Block buffer

(5% Nonfat milk in TBS-T)
5g milk in 100 ml TBST

10. NaCl 4 M

2 L: 467.5 g NaCl. Autoclavable.

11. NaOH 10 M

0.5 L: 200 g

12. NaAc 3 M

500 ml: add 204 g NaAc.3H₂O (FW 136), adjust pH by glacial acetic acid (~60 ml) to pH5.2. Autoclavable.

13. MgCl₂ 1M

500 ml: Add 101.65 g MgCl₂.6H₂O into 500 ml ddH₂O. Autoclavable.

14. CaCl₂ 1M

400 ml: Add 58.8 g CaCl₂.2H₂O (FW 147), filter for sterilization.
Dilute 10x to make 100 mM CaCl₂.

15. MgSO₄ 1M

500 ml: Add 123.3 g MgSO₄.7H₂O into 500 ml ddH₂O. Autoclavable.

16. ZnCl 0.5M

50 ml: 3.4 g ZnCl to 50 ml.

Stock in -20°C

1. IPTG (1 M)

1 g IPTG (FW 238.3) resolved in 4.2 ml (~4 ml) ddH₂O, filter through 0.22 µm filters, aliquot 1 ml in eppendorf. Store at -20°C.

2. DTT (1 M)

5 g DTT (FW 154.25) resolved in 32.5 ml (~30 ml) 10 mM NaAc (pH 5.2), filter through 0.22 µm filters, aliquot 1 ml in eppendorf. Store at -20°C.

3. X-gal (20mg/ml)

Add 5 ml (~4.8 ml) **DMSO** into 100 mg X-gal bottom (FW 408.24). Store at -20°C.

4. PMSF (100 mM, =17.4 mg/ml)

Resolve 1.74g PMSF (MW 174) in **isopropanol**, total 100 ml. Aliquot and store at -20°C or R.T..

5. Carbencillin or Ampicillin (50 mg/ml) in water. 1000x

2.5 g 50 ml.

6. Kanamycin (10 mg/ml) in water. 200x

0.5 g 50 ml.

7. Chloramphenicol (34 mg/ml) in ethanol. 200x

1.7 g/ 50 ml.

8. lysozyme 50 mg/ml, 1000x.

2.5 g/ 50 ml.

9. TSA (MW 303):

Add 1.32 ml Ethanol into each vial (1 mg?) to make the TSA stock 2.5 mM, 5000x. Final concentration of TSA in the cell culture is 0.5 µM (~150 ng/ml).

Solutions.

1. Bacteria lysis buffer (GST pull-down binding buffer)

(50 mM Tris 7.5, 150 mM NaCl, 0.05% NP-40.)

1L
50 ml 1M Tris HCl 7.5;
37.5 ml 4 M NaCl;
5 ml 10% NP-40.
ddH₂O to 1L.

1.2. GST pull-down binding buffer (1 M)

(50 mM Tris 7.5, 1 M NaCl, 1% NP-40.)

500 ml
25 ml 1M Tris HCl 7.5;
125 ml 4 M NaCl;
50 ml 10% NP-40.
ddH₂O to 500ml.

2. RIPA Buffer

(50 mM TrisHCl pH7.4, 150 mM NaCl, 2 mM EDTA, 1% NP-40, 0.1% SDS)

1L
50 ml 1 M Tris 7.4,
37.5 ml 4 M NaCl,
4 ml 0.5 M EDTA,
10 ml NP-40.
10 ml 10% SDS.

3. Cell Lysis Buffer (Flag-IP buffer)

(50 mM TrisHCl pH7.4, 250 mM NaCl, 0.5% Triton X100, 10% glycerol, 1 mM DTT, PMSF, PI (Roche))

1L
50 ml 1 M Tris 7.4,
62.5 ml 4 M NaCl,
5 ml Triton X-100,
1 ml 1 M DTT,
100 ml glycerol.

4. H-Lysis solution:

(0.25M sucrose (MW=342), 3 mM CaCl₂, 1 mM Tris pH8.0, 0.5% NP-40)

500 ml
43 g sucrose
1.5 ml 1 M CaCl₂
0.5 ml 1 M Tris pH8.0
25 ml 10% NP-40
add ddH₂O to 100 ml

Filter sterilize, store at 4°C.

5. H-Wash solution:

(300 mM NaCl, 5 mM MgCl₂, 5 mM DTT, 0.5% NP-40)

500 ml
37.5 ml 4 M NaCl
2.5 ml MgCl₂
2.5 ml DTT
25 ml 10% NP-40

6. H-Extraction solution (Histones Extraction):

(0.5 M HCl, 10% glycerol, 0.1 M 2-mercaptoethylamine-HCl (MW: 113.6)).

50 ml
2.25 ml HCl (11.2 M)
10 ml 50% glycerol.
0.55 g 2-mercaptoethylamine-HCl (Sigma name: cysteamine hydrochloride)

Recipe of making SDS-PAGE

SDS-PAGE		12% resolve gel	10% resolve gel	4% stacking gel 10 ml
20 ml (For 4x1mm plate)	4x buffer	10	5	2.5
	40% Acr-Bis	6	5	1
	ddH ₂ O	4	10	6.5
	10% APS (ul)	100	100	100
	TEMED (ul)	20	20	20
20 ml (For 4x1mm plate)	4x buffer	5	5	2.5
	30% Acr-Bis	8	6.6	1.3
	ddH ₂ O	7	8.4	6.2
	10% APS (ul)	100	100	100
	TEMED (ul)	20	20	20
	Voltage	Time	Buffer	Volume
Electrophoresis	150V/200V (15/30 mA)	1 h	Tris/Glycine/SDS	300 ml tank
Semi-Dry Transfer	5V (40-60 mA/gel)	1-2 h	20% Methanol 1x SDS running buffer	Make 500 ml for 4 gels
Agarose	100V	30min-1h	10 µl EB to 100 ml agarose/TAE	350 ml tank