

### Mononucleosome binding assay

**2x Binding buffer:** 40 mM HEPES 7.9, 160 mM KCl, 0.2 mM ZnCl<sub>2</sub>, 0.2% EDTA, 20% glycerol, 0.2% NP40, 1mM DTT, 1 mM PMSF.

Total vol: 20  $\mu$ l.

Mononucleosome: ~ 1-1.5  $\mu$ g. (4-5  $\mu$ l of 0.3 mg/ml) per binding.

Proteins: variant. Depend on binding affinity to nucleosome.

Core mix for 10 rx: 40  $\mu$ l mononucleosome + 40  $\mu$ l 2x binding buffer, aliquot 10  $\mu$ l each.

Add proteins: GST-ING2PHD: 0, 5, 10, 20, 50  $\mu$ g (1  $\mu$ l, 2 $\mu$ l, 4  $\mu$ l and 10  $\mu$ l of 5 mg/ml).

GST control 0, 20, 50  $\mu$ g (5 $\mu$ l and 12.5  $\mu$ l of 4 mg/ml). (diluted with binding buffer).

Binding condition: leave at 30°C 30 min.

Add 5  $\mu$ l 6x DNA sample buffer, load onto 5% acryamide/TBE gel, ran with **1xTBE buffer** under 120 V 60-70 min.

Stained in EB/TBE (5  $\mu$ l EB in 50 ml TBE) for 1 h.