

Chromatin Isolation by small-scale biochemical Fractionation

Mendez, J., and B. Stillman. 2000. Mol. Cell. Biol. **20**:8602-8612.

Buffer A

10 mM HEPES, pH 7.9
10 mM KCl
1.5 mM MgCl₂
0.34 M Sucrose
10 % Glycerol
1 mM DTT
Protease inhibitor cocktail

Buffer B

3 mM EDTA
0.2 mM EGTA
1 mM DTT
Protease inhibitor cocktail

1. Harvest $1 \times 10^7 - 2 \times 10^7$ cells by using a cell scraper; spin down at 1000 rpm for 2 min, discard supernatant
2. Wash cell pellet with PBS, spin down at 1000 rpm for 2 min.
3. Repeat second step.
4. Resuspend cell pellet in 200 μ l of Buffer A.
5. Add Triton X-100 to a final concentration of 0.1%.
6. Incubate cells on ice for 8 min.
7. Centrifuge at 1,300 x g, 4 °C, for 5 min; separate supernatant = fraction **S1** from pellet (nuclei) = fraction **P1**.
8. Clarify **S1** by high-speed centrifugation at 20,000 x g, 4 °C, for 5 min; collect supernatant= fraction **S2** (discard P2)
9. Wash **P1** once with Buffer A and lyse it for 30 min in Buffer B (100 μ l).
10. Centrifuge at 1,700 x g, 4 °C, for 5 min; separate supernatant = fraction **S3** from pellet (chromatin) = fraction **P3**.
11. Wash **P3** once with Buffer B and resuspend it either in SDS sample buffer (then boil for 10 min at 70 °C and analyse chromatin associated proteins by SDS PAGE/Western Blot) or in nuclease digestion buffer.